# Sulfide



# test paper

Sulfide test paper is an equal alternative to lead acetate paper, but with an important advantage, that it does not contain any hazardous materials, subject to declaration.

### Colour reaction:

This white test paper changes, in the presence of sulfide ions and hydrogen sulfide, to brown-black, depending on the concentrations.

#### Instructions for use:

Place a drop of the sample on the test paper, which will then change colour to brown-black. Smaller quantities of sulfide now appear as a brown ring.

For the determination of hydrogen sulfide, first wet the test paper with distilled water, and hold it up in the air to be tested. A brown colouration, on that part of the test paper where water was applied, shows the presence of hydrogen sulfide.

# Limit of sensitivity:

The presence of 5 mg/l sulfide (S<sup>2-</sup>) produces an easily recognizable brown ring.

MACHEREY-NAGEL GmbH & Co. KG Neumann-Neander-Str. 6-8 · D-52355 Dueren (Germany) Phone +49 2421 9 69-0 · Fax +49 2421 9 69-199 e-mail: sales-de@mn-net.com · www.mn-net.com